

ECO 5 NEW

Ultra-Portable, Ultra-Affordable, Color doppler



CHISON
Value Beyond Imaging

Eco-friendly

- Less radiation
- Lower power consumption
- Longer life: LED

Ergonomic


- Rotatable LED monitor (0~30°): tailor to different viewing angles
- Unique LED: reduce eye fatigue
- Keyboard design: for easy and quick access
- Full screen display mode



Excellent

- Full digital beam-former
- Color doppler technology
- Double frame rate technology
- Color priority
- Color filter
- Applications
- Abdomen, OB&GYN, Urology, Vascular,
- Small Parts, Pediatrics, MSK, Basic Cardiology
- Display modes: B, B/B, 4B, M, B/M, CFM, PW, Trapezoidal
- Advanced technologies: i-Image, THI, SRA, Compound imaging
- 16G memory card
- 320G harddisk (option)
- streamlined workflow





• **Streamlined Workflow**


The multi-purpose user presets, comprehensive measurement & report system, built-in EasyView™ image archive system, quick image storage / retrieve / transfer, one-button direct print, make the complete workflow better than what you can dream of.



Various Measurements



Detailed Reports



EasyView™ Image Archive System

• **ONE KEY**



• **Support PC Printer (HP1102&M251n)**



2.5MHz-5.0MHz Convex
C3-A



5.3MHz-10.0MHz Linear
L7M-A



5.3MHz-11.0MHz Linear
L7S-A



2.5MHz-5.0MHz Micro-Convex
MC3-A



4.5MHz-8.0MHz Transvaginal
V6-A



4.5MHz-8.0MHz Micro-Convex
MC6-A



5.0MHz-10.0MHz Transrectal
R7-A

CHISON MEDICAL IMAGING CO., LTD.

No. 9 Xin Hui Huan Road, New District, Wuxi, Jiang Su Province, China214028
TEL: 0086-510-85310593 / 85310937 FAX: 0086-510-85310726 EMAIL: export@chison.com.cn



We reserve the right to make changes to this catalogue without prior notice.
Please contact us or your local dealer for the latest information.
* Please refer to the user manual for details

ECO5-141028