

CHISON
Value Beyond Imaging

SonoTouch 30

NEW



Get in **Touch** with
The Future of **Ultrasound**

Approved by FDA & CE

Ready to GO

GoLED™

- The entire button-free screen is touch-driven
- Waterproof front panel is easy to sanitize

GoMeasure™

- Make measurements with your fingers only
- Measurement accuracy is under 1 mm

GoHub™

- Quick copy, paste and transfer of images & reports
- Easy import and export of patient information or case studies
- USB, DICOM 3.0,(Option)

GoPower™

- Quick boot in 30 seconds
- Car charger available
- Instant shutdown in 2 seconds

GoTune™



- Adjust image settings like TGC, focus and depth with intuitive touch-driven instructions

GoApp™



- Select your application with visual icons

GoFlex™

- Adjustable tilt angle
- Make your measurement easier
- Also functions as a convertible handle for carrying

GoClip™

- Quick save and recall
- Plays CINE clips forward or backward



PROFESSIONAL Clinical Application UI

The dedicated icon-based presets for major MSK / Anaesthesia/ Vascular positions will help doctors acquire optimal images in an easier, quicker, and more intuitive way.


It uses the latest Super Tissue technology on the 11MHz Linear probe to enhance the near field image within 3 cm depth which is the focal region of most Point of care applications.



| | | | |
|--------------------|--|---|---------------------------------------|
| <p>MSK</p> | <p>Subscapularis Tendon, B Mode</p> | <p>Patella Tendon Tibian Insertion Part, B Mode</p> | <p>Palmaris Longus Tendon, B Mode</p> |
| <p>Vascular</p> | <p>Brachial Artery, B Mode</p> | <p>Carotid, PW Mode(Auto trace)</p> | <p>Central Venous Access, B Mode</p> |
| <p>Anaesthesia</p> | <p>Popliteal Sciatic Nerve, B Mode</p> | <p>Median Nerve</p> | <p>Anchillies Tendon, B Mode</p> |

Super Needle (Option)

Super Needle (Option) is our needle enhancement technology supporting angle adjustment up to +/- 30 degrees. It enhances the needle's appearance for a more confident approach to anatomical targets.



Without Super Needle
Needle is Not clear

With Super Needle
Needle is Better Visualized