

Perfect for small animal examinations


- Exceptional image quality
- 10.4" touch screen monitor
- Compatible with different probes
- Mains or battery power supply
- Easy to disinfect



DRAMIŃSKI 4vet

DRAMIŃSKI 4vet

PORTABLE VETERINARY ULTRASOUND SCANNER WITH ELECTRONIC PROBES



LCD touch screen

Digital imaging technology. (B, B+B, B+M) for thorough examinations of small and large animals in a mixed veterinary practice.

- 10.4" LCD screen with LED backlight
- Clear and high contrast image
- Gamma, Gain, Dynamic Focus

Saving to hard drive

Windows operating system. Easy image and Cine loops saving and exporting.

- Internal hard drive
- Internal PC with Windows operating system
- Communication ports: USB 2.0, LAN, VGA

Perfect for all vet clinics

Mains or battery power supply.



- Battery: Li-ion 6.2 Ah
- Operation time on a single battery: approx. 2 hours
- Optional stand for stationary use

Duralumin housing

DRAMIŃSKI 4vet casing allows easy disinfection of the device. Full range of different electronic probes.



- Touch panel and user friendly menu in different languages
- Broadband probes: microconvex 4-9 MHz, linear 40 mm 6-14 MHz, convex 2-8 MHz